

AIR LECTURE 9 & 10

Sizing of Particles

When dealt with powders \rightarrow you sieve the particles

You take a vessels with holes in it (of a certain size).

You put the particles and start shaking it.

Using a series of sieves is a common practice done by chemical and civil engineers known as sieve analysis.

Assumption \rightarrow Particle has a certain gravitational settling rate.
ie there will be enough velocity due to gravity.

This however isn't the case for particles ^{smaller} ~~lower~~ than 100 micrometer.

Buoyancy/ Drag makes sieving impractical and of no use in particles of this small size.


How do you size/count particles in nm size range?

\rightarrow Inertial Impactors (MOUDI)

Micro orifice Uniform
Deposition Impactor.

Diameter of orifice becomes smaller as
you go further down.

Using fluid mechanics to separate smaller &
large particles.



Working of MOUDI

- 1) There are 3 sizes of particles. They enter from atmosphere.
- 2) There is a pump to create vacuum (suck atmosphere out)
- 3) They come in with certain velocity and hence an inertia.
- 4) The gas molecules (the air) adapt to the change caused at the impaction stage (obstruction).
- 5) The gas particles are forced to turn.
- 6) The smaller particles turn away, the larger ones are unable to turn and get stuck to the vaseline on impaction stage.
- 7) This process continues for every size of particles.
- 8) Now you can get histogram of data collected.

In industry, only 2 stage process, $pm-10$ & $pm-2.5$.

Particulate matter

→ Used for particles ($> 56 nm$) and is mass based.

Optical Particle Counter (OPC's)

- Number based ($> 100 nm$)

You count the particles one by one.

→ Light scattering / Extinction


Limited to $0.09 \mu m$.

$$0.4 > \lambda_{opt} > 0.7 \mu m$$

→ Aerosol → single particle
Water → Multiple particles.

Cloud of particles
↓
turbidity.

↓
We measure
turbidity in
water.



- Larger the particles → longer the signal
- Every time a particle passes, there is a beep.
- Each beep is a count = size of particles.

Electrical Mobility

- Sizing ($> 3 \text{ nm}$)
- Counting (Condensation Nuclei Counters)
- Use an electrical field to charge particles and then you use an electric field to separate them

Working (Model 3081 Long DMA) → fig in next page.


- 1) You bring in the particle with a charge and you introduce them into the shell.
- 2) There is an electric field between the outer surface and the brown (high-voltage) rod.
- 3) The ~~charge~~ charged particle moves under the influence of electric field.
- 4) Depending on the size, there will be a charge and depending on the charge, there will be an electrical mobility.
- 5) The dotted lines (can be seen faintly) are the trajectories that indicate the different size of particles that hit the central rod.
- 6) Below the rod, there is a small opening, which allows only one particular trajectory / size to go through.
- 7) With a particular electrical field (tuned electrical field), a particular size can go through the opening

Polydisperse particles → Monodisperse particles.


To change size of monodisperse particle → change voltage diff between central rod & shell

Electrical Mobility Analysers

- Charge particles
- Electrical Field



SIZE




CPC

**OPC in Principle – Alcohol
used to grow size**

**MODEL 3081
LONG DMA**

Differential Mobility Analyser (For size separation using electrical mobility)



Source : TSI Inc.

Now to count the particles

↓
Condensation Particle Counter. (3 nm)

- 1) Pass the monodisperse particles through an area saturated with vapor/alcohol. These vapor/alcohol condenses on the small particle to create a larger particle which is optically detectable. ~~this is called scanning mobility particle sizer.~~

Scanning mobility Particle Sizer.

Combination of differential mobility analyzer (sizer) and a counter (CPC).

Long unit is has the electrical rod.

Why scanning? You scan the electric level, i.e. you choose electric field as per your demand.

Field Instruments

Standard sampler → High volume sampler.

(Govt used)



Take large volume of air & pass it through a filter which has a cyclone cutoff (removes particles $> 10 \mu m$).

So you are only collecting particles that are respirable on a filter paper. Weigh the filter paper before and after the collection of particles.

W_{after} - W_{ini}


W_{particles}

↓ Div by Volume

$\frac{W_p}{V_{\text{sample of air}}}$

Condensation Particle Counter (For growing the size separated particles and detecting them using light scattering)

MODEL 3022A CONDENSATION PARTICLE COUNTER



Scanning Mobility Particle Sizer



Field Instruments

- High Volume Samplers (~2 Lakhs)
 - Regulatory (RSPM, SO_x, NO_x)



- Mini Vol Samplers (5 l instead of 1000 l)
 - For places that are highly ~~poll~~ polluted.
 - Principle used: same as MOUDI

→ Dust Trak

You put an impactor at entry & measure particles less than the given size.

Measuring Gaseous Pollutants

• Wet ~~chemistry~~ chemistry

Collect the sample → Take a solvent reagent → Bubble the air for 8-24 hrs → Bring the air back to lab → Perform standard test protocols to measure oxygen/sulfur/nitrogen oxides.

• Electrochemical sensors (Real time)

→ Instant

→ Accurate (not as much as wet chemistry) but more feasible.

• Spectroscopy

Field Instruments

- High Volume Samplers (~2 Lakhs)
 - Regulatory (RSPM, SO_x, NO_x)
- Mini-vol Samplers (~5 Lakhs)
 - Regulatory (PM₁₀, PM_{2.5})



Field Instruments

- High Volume Samplers (~2 Lakhs)
 - Regulatory (RSPM, SO_x, NO_x)
- Mini-vol Samplers (~5 Lakhs)
 - Regulatory (PM₁₀, PM_{2.5})
- Dust Trak (~7 Lakhs)
 - Real time (1 minute resolution)
 - PM_x (x = 1, 2.5, 4, 10)



Measuring Gaseous Pollutants

- Wet Chemistry (8 hour averages)
- Electrochemical sensors (Real time)
- Spectroscopy (Real time)
- Standard Methods (SPCB's Analytical Labs)